

EASYADAPT EVOLV SIT-TO STAND FRAME

Modular design makes the EasyStand Evolv the most versatile standing frame available. The base unit is functional for many users, however when more support is needed a wide variety of positioning options can be added to accommodate the most involved users.

Medical Benefits of Standing

- Improve/maintain range of motion
- Decrease joint/muscle contractures
- Management of atrophy in the trunk and leg muscles
- Improve strength to trunk and lower extremities
- Decrease muscle spasms
- Improve/maintain bone integrity/skeletal development
- Lessen/manage the progression of scoliosis
- Manage pressure (ulcers) through changing positions
- Improve bowel function and regularity
- Aid in kidney and bladder functions
- Strengthen cardiovascular system and build endurance
- Improve circulation
- Reduce swelling in lower extremities



Standard Front
with options



Shadow Tray
with options

NAME	EASY STAND EVOLV SIT-TO-STAND
PRODUCT CODE	PNG50084
CONTROLS	MEDICAL HAND PUMP
CASTORS	6" SOLID
STD ACCESSORIES	SEAT – BACK AND TRAY UPCHARGE OPTIONS
LOAD CAPACITY	LOAD CAPACITY - 115KG
ADJUSTABILITY	SEAT DEPTH

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za