

## PHOENIX STAND-UP POWER



<b>NAME</b>	SATURN PHOENIX POWER
<b>PRODUCT CODE</b>	49312
<b>SIZE</b>	14", 16", 18"
<b>BATTERIES</b>	2XSEALED 36AH/12 VOLT
<b>CASTORS</b>	8X2"FOAM FILLED GREY TYRE
<b>REAR TYRES</b>	3.00-8 FOAM FILLED GREY TYRE
<b>CONTROLS</b>	DYNAMIC SHARK REMD11 DUAL FUNCTION SHARK REMOTE REMD11 POWER MODULE PMB11
<b>LOAD CAPACITY</b>	LOAD CAPACITY – 110KG GROSS WEIGHT – 90KG
<b>ADJUSTABILITY</b>	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT
<b>STD ACCESSORIES</b>	SOLID FIXED CLOTHES GUARDS PERSPEX
<b>PACKING</b>	103X72X72CM

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za