

Power Wheelchair Outdoor/Indoor

Invacare Stream Power
CTM Rehab Power (Standard, with Tilt, Bed Type)
TDX SP2NB Power
Invacare Kite Power
Storm 4
CTM Bariatric Power

INVACARE STREAM POWER



Maintain independence

The **Stream**, a new power wheelchair from Invacare is a popular model suitable for with reduced mobility. The Invacare **Stream** is ideal for both indoor and outdoor use. To ensure a comfortable sitting position and a relaxing ride, the seat can be adjusted to suit customers unique requirements. A dedicated lighting kit is available as an option, thus extending the possibilities of use for this model.

PRODUCT CODE: 1533501



Light-kit
Available as an option - includes bright headlights with indicators.



Rear wheels
14" pneumatic rear wheels with studded profile.



Swing in/out footrest
The leg rests swing in and out for easy transfer and transportation.

Technical data

Stream	44 - 49 cm ¹	41, 46, 51 cm	42,5 - 50 cm ²	48 cm ²	29 - 36 cm ²	30 - 43 cm ²	up to 22°
Stream	-10° - 30°	64 cm ³	100 cm	103 cm	105 kg	130 kg	50 Ah
Stream	166 cm	5 cm	11%	2 x 180 W	25 km ⁴	Frame colour Gloss aluminium	Upholstery Black nylon
Stream	6 or 9 km/h	<ol style="list-style-type: none"> 1/ seat cushion 43 cm 2/ measured without seat cushion / upholstery 3/ total width is depending on seat adjustment 4/ according to ISO 7176-4:2008 					

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CTM REHAB POWER (STANDARD, WITH TILT, BED TYPE)



Lights & cushion not included



Features

- * Maximum Speed: 8kmph/5mph
- * Weight Capacity:
7200: 7200T: 162kg / 357lbs
7200B: 136kg / 300lbs
- * Suspension System
- * Standard 75ah Batteries and Easy Accessible Design
- * Manual Tilting and Adjustable Back Angle
- * Height Adjustable / Flip Up Armrests

NAME	CTM REHAB POWER
PRODUCT CODE	HS7200T
SEAT WIDTH	16", 18", 20"
BACKREST	ANGLE ADJUSTABLE
CASTOR	9" PNEUMATIC
REAR WHEEL	14" PNEUMATIC
BATTERIES	SEALED 50 A/H
ARMS	FLIP BACK HEIGHT & WIDTH ADJUSTABLE
FOOTREST	REMOVABLE
FRAME	FIXED MILD STEEL
OPTION	TILT IN SPACE


Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

TDX SP2^{NB} POWER

The remarkable *Invacare TDX SP2^{NB}*, with its **compact narrow base** boasts **durability** and **ultimate performance**. Great manoeuvrability indoors and out, with a wide choice of seating for **optimum comfort and support** makes the *TDX SP2^{NB}* the powerchair for individuals who need full functionality, but do not want to compromise on style. Key features include:

- Enhanced stability
- Intuitive driving experience
- Simply Smart LiNX control system



Technical data

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle
<i>Modulite™</i> 1-piece seat plate	440 - 480 mm ¹	410 - 510 mm	450 mm ⁴	480 mm ⁷	275 - 340 mm	290 - 460 mm ⁸	0 - 8°
<i>Modulite™</i> telescopic seat	380 - 530 mm ²	410 - 510 mm 510 - 560 mm ³	450 / 475 mm ^{4,5} 475 - 775 mm ⁶	480, 540, 560 mm ⁷	245 - 310 mm	290 - 460 mm ⁸	0 - 20° / 0 - 30° ⁵
	Backrest angle	Total height	Driving unit width	Total length inc. Footrests	Total length without footrests	Max. User weight	Battery capacity
<i>Modulite™</i> 1-piece seat plate	90° - 120°	1010 mm ¹¹	610 mm ¹⁰	1160 mm ⁹	980 mm	150 kg	50 Ah
<i>Modulite™</i> telescopic seat	90° - 120°	1010 - 1170 mm ¹¹					
	Turning diameter	Max. Climbable obstacle height	Max. Safe slope	Motor capacity	Speed	Driving range	Total weight
	1120 mm	75 mm	9° / 16 %	2 x 340 W	up to 26 km ¹²	6 and 10 km/h	138 kg ¹¹



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE KITE POWER



The *Invacare Kite* is a powerchair that combines **performance, driving comfort, compactness and personalisation**. Designed for **active users**, this remarkable chair has a **compact design** making it ideal for outdoor activities and comes with a **wide range of colour options**. Its **Dual Swing Technology (D.S.T.)[®]**, provides a unique suspension system to guarantee a **superior driving experience with improved comfort and traction**. The optional **smooth-drive front castors** are set to further improve your comfort when out and about, especially over uneven terrain.

- Patented Dual Swing Technology (D.S.T.)[®]
- Superior drive abilities
- Configured for customisation
- LiNX Smart Technology



PRODUCT CODE: KITEAA2

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Seat angle	Total weight
One piece seat plate	440 - 480 mm ¹	410 - 510 mm	440, 460, 480, 500 mm ³	480 mm ³	275 - 340 / 325 - 390 mm	0 - 20°	125 kg ⁹
Modulite [™] telescopic seat	380 - 530 mm ²	410 - 510 mm	460, 480, 500 mm ³	480, 540 or 560 - 660 mm ^{3, 6}	245 - 310 / 295 - 360 mm	0 - 20° / 0 - 30° ⁵	126 kg ⁹
Max Seat	560 - 600 mm	510 - 580 mm	460, 480, 500 mm ³	630 mm ^{3, 14}	290 - 360 mm	0 - 12° / 0 - 20° / 0 - 28° ⁵	127 kg ⁹
	Length of footrests	Backrest angle	Driving unit width	Total height	Total length inc. footrests	Total length without footrests	Max. user weight
<i>Invacare Kite</i>	290 - 470 mm ^{3, 4}	90° - 120°	595 mm ⁷	max. 1170 mm ⁸	max. 1140 mm ⁸	max. 870 mm ⁸	136 / 160 kg ¹⁰
	Battery capacity	Turning diameter	Max. climbable obstacle height	Max. safe slope	Motor capacity	Driving range	Speed
<i>Invacare Kite</i>	60 Ah / 73 Ah	ca. 1600 mm	60 mm ¹¹	10° / 18% ¹²	2 x 250 W / 2 x 340 W ¹⁰	up to 39 km ¹³	6 / 10 km/h

28° Tilt Mode

The 28° Tilt mode with centre-of-gravityshift is for those who require a greater repositioning tilt angle for optimum posture



Standard Seat Tilt

Seat angle can be set by a spindle or adjusted electronically for greater stability and improved downhill comfort.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

STORM 4

The *Invacare Storm⁴* powerchair was created to satisfy complex requirements regarding **configurability, adaptability and functionality**.

The *Storm⁴* has impressive outdoor and indoor performance providing a safe and stable drive experience without compromising on manoeuvrability. Thanks to its contemporary design, it also ensures the powerchair is always elegant and stylish. Key features include:

- **Highly configurable**
- **Easy handling**
- **Aesthetic design**
- **LiNX Smart Technology**



Technical data

For more comprehensive clients information about this product, including the product's user manual, please visit www.invacare.co.uk

	Seat width	Seat depth	Backrest height	Total height	Backrest angle	Armrest height	Total weight
<i>Modulite™</i>	380 - 530 mm	410 - 510 mm	480, 540, 560 mm ^{1,10}	1020 - 1190 mm	90 - 120°	275 - 340/ 325 - 390 mm ¹¹ 245 - 310/ 295 - 360 mm ¹²	174 kg ³
<i>Recaro</i>	490 - 530 (360) mm	460 - 510 mm	770 - 830 mm ¹	1300 mm	90 - 135°	250 - 340 / 290 - 380 mm	174 kg ³
	Seat height	Length of footrests	Seat angle	Driving unit width	Total length inc. footrests	Total length without footrests	Max. user weight
<i>Invacare Storm⁴</i>	400 - 650 mm ⁴ 450 / 480 mm without lifter	290 - 460 mm ⁵	0° - 25° ⁹	630 mm (wheelbase)	1190 mm ⁷	910 mm	150 kg
	Battery capacity	Turning diameter	Max. climbable obstacle height	Max. safe slope	Motor power	Speed	Driving range
<i>Invacare Storm⁴</i>	73.5 Ah	1770 mm	60 / 100 mm	20% (11.3°) 15% (8°) with lifter	340 W	6 / 10 km/h	up to 36 km ⁸



Seat tilt

Manual settable or powered tilt with maximum 25° for optimum comfort and posture.



Telescopic seat frame

This frame is the heart of the *Modulite* seat design, supporting a number of positioning options.



Driving unit design

Attractive and modern design with easy to change, colourful shroud covers.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CTM BARIATRIC POWER



Solid Tyre Option

Features:

- Maximum Speed: 8kmph/5mph
- Capacity: 181 kg / 400lbs
- 500mm / 20" reclining captains chair with seat belt
- Swing-away joystick controller
- ZERO effort single free wheel lever

PRODUCT CODE: HS5600



HS-5600			
Type	Power Chairs	Curb Climbing	50 mm / 2"
Model Number	HS-5600	Turning Radius	790 mm / 31.1"
Overall Length	1135 mm / 44.7"	Supension	N/A
Overall Width	635 mm / 25"	Brake	Electro-Mechanical
Overall Height	1185 mm / 46.7"	Seat Type	Captain Seat, Adjustable Seat Height, & Seat Belt
Wheels: Front	230 mm / 9"	Seat Width	500 mm / 20"
Wheels: Rear	360 mm / 14"	Motor Size	420W, 4600 r.p.m
Weight w/ Batteries	102 kg / 224 lbs	Battery Size	(2) 12V. 50Ah
Max. Speed	8 kmph / 5 mph	Battery Weight	31.3 kg / 69 lbs
Weight Capacity	181 kg / 400 lbs	Travel Range	24 km / 15 Miles
Ground Clearance	85 mm / 3.4"	Battery Charger	5A Off Board
Grade Climbable	8 degree	Electronics	N/A
Color	Burgundy, Silver, Blue		

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za