

# CHASSWHEEL FOURX

## GOING PLACES

Chasswheel produce top of the range all terrain power chairs for individuals who choose an active lifestyle. Its innovative design features 4-wheel suspension and flexible chassis that makes it ideal for challenging terrain. With Chasswheel, going off the beaten track provides an exciting, all round driving experience.



## TECHNICAL SPECIFICATIONS

Width:	695 mm
Height (normal seating position w/out headrest):	700-850 mm
Length:	1250 mm (footrest down)
Weight with battery:	140 Kg
	Four X DL Lift 170 Kg
	Four X DL SSS 170 Kg
Maximum speed:	10 Km/h
Drive distance:	~40km
Batteries:	2 x 73Ah, gel
Motors:	4 x 250W
Max user weight:	125 Kg (160Kg with HD kit)
Seat width:	400-480 mm
Seat depth:	390-510 mm
Backrest height:	500-600 mm
Ground clearance:	6" tyres 100mm
	8" tyres 120 mm
Slope climbing capacity:	20°
Maximum obstacle climbing ability:	150 mm
Seating height without cushions;	6" tyres - Four X DL 440 mm
	Four X Lift - 530 mm
	Four X SSS - 530 mm
	8" tyres - Four X DL 465 mm
	Four X Lift - 555 mm
	Four X SSS - 555 mm
Gravity point adjustment:	250 mm
Tilt angle:	30°
Footrest angle:	80°
Backrest angle:	84°-125°
Turning radius:	1.15m
Electronics:	Four X, Four X DL, Four X Lift
	Dynamic DX2
	Four X SSS - R-Net



# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za