

# CTM BARIATRIC SCOOTER



**Solid Tyre Option**



PRODUCT CODE: HS890

HS-890			
<b>Type</b>	Mobility Scooters	<b>Curb Climbing</b>	70 mm / 2.8"
<b>Model Number</b>	HS-890	<b>Turning Radius</b>	1730 mm / 68.1"
<b>Overall Length</b>	1530 mm / 60.2"	<b>Supension</b>	Front & Rear
<b>Overall Width</b>	685 mm / 27"	<b>Brake</b>	Hand Brake & Electro-Mechanical
<b>Overall Height</b>	1245 mm / 49"	<b>Seat Type</b>	Swivel Captain Seat With Seat Sliding Mechanism
<b>Wheels: Front</b>	320 mm / 13"	<b>Seat Width</b>	480 mm / 19"
<b>Wheels: Rear</b>	320 mm / 13"	<b>Motor Size</b>	700W 5400 r.p.m
<b>Weight w/ Batteries</b>	114 kg / 251 lbs	<b>Battery Size</b>	(2) 12V. 50Ah
<b>Max. Speed</b>	15 kmph / 9.3 mph	<b>Battery Weight</b>	31.3 kg / 69 lbs
<b>Weight Capacity</b>	227 kg / 500 lbs	<b>Travel Range</b>	37.5 km / 23.3 Miles
<b>Ground Clearance</b>	85 mm / 3.3"	<b>Battery Charger</b>	5A Off Board
<b>Grade Climbable</b>	10 degree	<b>Electronics</b>	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
<b>Color</b>	Burgundy, Silver, Blue, Turquoise		

\* Subject to change without notice.

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za