

INT ROADSTER BUGGY



Yk8001 which seat width is 280mm-320mm
 YK8002 which seat width is 320mm-400mm

- Aluminium frame, light weight and luxury design
- Foldable chair frame for easy transportation and storage
- Comfortable soft seat, back over and adjustable soft headrest
- Height adjustable & foldable Pushbar
- Removable solid base seat system
- Seat Back height adjustable and seat base depth adjustable
- Reclining Backrest
- Adjustable angle & height footrest
- Mesh bag under for small belongings

Seat width 280mm-320mm YK8001
 Seat width 320mm-400mm YK8002
 Color: Orange & grey
 Blue & grey



FOLDED & SEAT REMOVED

FOLDED:	Specification 82cm*66cm*45cm	Yk8001	20kg
FOLDED:	Specification 82cm*74cm*45cm	Yk8002	22kg