

# INT SPRY ACTIVE WHEELCHAIR



Code: YKACTIVE

Aluminium Powder Coated Frame  
Front Castors 5" solid  
QR L/Weight Rear Spoke Wheel with PU tires  
Adjustable COG  
Detachable & Swing away armrest  
Single Piece footplate height adjustable  
Seat and Backrest Black Nylon  
TAB back Upholstery  
Integrated height adjustable push handle  
Flip Down & angle adjustable Backrest  
Undermount Brake system  
OPTION: Anti-tipper  
Sizes: 12", 14", 16", 17", 18", 20"  
Seat Depth: 16" & 18" (seat width dependent)

Carton size: 68x62x60 CM  
GW/NW: 13/12 KG

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za