

TIAN FEEDING & POSTURE CHAIRS

PRODUCT CODE: TC401 & TC403



Code: TC401

SPECIFICATION

1. Dimensions(L*W*H)mm:630*430*(840~950).
2. Seat size(L*W*H)mm:330*(300~360)*360.
3. Load capacity: ≥ 70 kg
4. Net weight: 16kg
5. Adopting carbon dioxide welding, electrostatic spraying, anti-aging and good anti-rust.
6. All pad are a foam forming.
7. Distance can be adjust of pillow height: 0-140mm
8. Gear pad adjustment: 0-100mm
9. Pedal height adjustment:0-80mm.
10. Plastic moulded tray, quick & simple dismount for easy of use
11. The distance of front-back wheels: diameter ≥ 50 mm, Adjustable.

SPECIFICATION

- Table to seat cushion distance 210mm
- Head pad lift adjustment range 0~200mm
- Crotch pad before and after the adjustment range 0~100mm
- Seat to the maximum distance of the pedal 200mm
- Load Capacity 40kg
- Overall dimensions (L×W×H)680×560×800



Code: TC403