

# BUDDY BUGGY URBAN



Item Code Yk8003

Frame Powder coating  
Front wheel 8"×1.25" ,PU tires  
Rear wheel 12" PU tires  
Handgrip PVC  
Seat and backrest Nylon and woven cloth, Depth and width fixed  
Front fork Aluminum  
Footplate Plastic  
Anti-tipper tube  
Quick release for rear wheel  
Footrest height adjustable  
Breathable headrest and backrest  
High-quality sunshade"  
Seat Width 14"

Carton size: 1300\*420\*360mm  
GW/NW: 18.5/16.0KG

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za