



Code: Y205

Finding a folding frame aluminium Power wheelchair with large front castors and rear tyres has been a time old problem for wheelchair users particularly those in more rural areas. The solution has been neatly packaged and presented to the local market with the Bold looking Glide JUMPER Power wheelchair. Although not a full off-road the easily serviceable powerful motors and PG Electronics provide the user with assurance on quality and reliability allowing users to tackle uneven terrain in a power wheelchair that can be easily transported.



TECHNICAL SPECIFICATIONS

Product Code	Y205
Material	Aluminium
Sizes	Available in 16", 18", 19"
Wheels	Front Castor 9" foam filled Rear wheel 14" Foam filled
Battery	Lead acid battery 12V 34Amp
Motor	320watt*2
Slope	8 degrees
Armrest	Sloped Height Adjustable Arm
Controller	PG VSI 50A (UK)
Load Capacity	120kg
Colours	Red & Silver
Footrest	Swing Away Footrest
Seat	Detachable quality seat and back cushion, easy to clean
Special:	Advanced height adjustable anti-tip wheels. Highly reliable industry-proven electronics by PG Drives Technology
Max speed	5-8km
Travel mileage	20-25km
Weight of Wheelchair	40kg (excl battery - 20kg)

Overall Length (cm)	107
Overall Width (cm)	65
Overall Height (cm)	97
Overall Width Folded	36
Seat Width (cm)	43-48
Seat Depth (cm)	44
Back Height (cm)	44
Floor to seat Height	53
Weight Capacity (kg)	135
Front Caster Spec.	9"
Rear Wheel Spec.	14"
G./N. Weight (kg)	81.0/62.0
Carton Size (cm)	89x42x104



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za