

Seat Width Choices	18" (46 cm)
Seat Depth for 14" seat width	18" (46 cm)
Seat Depth for 18" seat width	20" (51 cm)
Seat Height	Front : 24.8"(63cm) Rear : 23.6" (60cm)
Backrest Height	62 cm (24.4")
Total Weight 14", 16",18"	284.5 lb/ 291 lb/ 298 lb (129 kg/ 132 kg/ 135 kg)
Open Size 14"(L x W x H)	45" x 25.5" x 48"(114 cm x 65 cm x 122 cm)
Open Size 16" (L x W x H)	45" x 28" x 48"(114 cm x 70 cm x 122 cm)
Open size 18" (L x W x H)	45" x 30" x 48(114 cm x 75 cm x 122 cm)
Footrest to Seat Length (Adjustable)	15.7" ~21" (40 cm ~ 54 cm)
Motor	DC24V, 320W
Battery	12V36Ah X 2 pcs

Packing Size 14" & 16" (L x W x H)	43" x 26" x 32" (110 cm x 67 cm x 82 cm)
Packing Size 18" (L x W x H)	43" x 28" x 32" (110 cm x 71 cm x 82 cm)
Front Wheel	9.5"
Mid Wheel	13"
Rear Wheel	8"
Slope Limitation	10°
Curb Limitation	2" (5 cm)
Ground clearance	2" (5 cm)

PHOENIX DRACO – POWER WHEELCHAIR

