

PIONEER POWER NG

Y201NG



Frame:	Reinforced Steel
Controller:	Yanteon Controller
Motor:	200W*2, brush motor with electromagnetic brake
Battery:	Lead-acid 2 x 12V 32Ah
Charger:	HP 24V, 3A, Taiwan
Mas speed:	7 km/h
Travel Mileage:	20km
Front wheel:	8" PU tire
Rear wheel:	12.5" PU tire, aluminium alloy
Armrest:	Hinged and height adjustable Armrest
Footrest:	Detachable, heel strap
Backrest:	TAB
Seat Width:	16", 18", 20"
Carton size:	87x38x81cm
G.W./N.W.:	63kg/55.3kg

		Steel	Aluminum
Overall Length (cm)		108	108
Overall Width (cm)		63	63
Overall Height (cm)		97	97
Overall Width Folded		34	34
Seat Width (cm)		43-48	43-48
Seat Depth (cm)		43	43
Back Height (cm)		45	45
Floor to seat Height		52	52
Weight Capacity (kg)		135	135
Front Caster Spec.		8"	8"
Rear Wheel Spec.		12.5"	12.5"
G./N. Weight (kg)		63.0/55.3	58.0/50.3
Carton Size (cm)		87x38x81	87x38x81